

# **PBGL CNV-seq Analysis v1.0**

## **A Laboratory Manual**

**Anibal E. Morales-Zambrana**

**Plant Breeding and Genetics Laboratory  
FAO/IAEA Joint Division  
Seibersdorf, Austria**

Created: March, 2021  
Last updated: 04 June 2021

**Please note:** *This is not an official IAEA publication but is made available as working material. The material has not undergone an official review by the IAEA. The views expressed do not necessarily reflect those of the International Atomic Energy Agency or its Member States and remain the responsibility of the contributors. The use of particular designations of countries or territories does not imply any judgement by the publisher, the IAEA, as to the legal status of such countries or territories, of their authorities and institutions or of the delimitation of their boundaries. The mention of names of specific companies or products (whether or not indicated as registered) does not imply any intention to infringe proprietary rights, nor should it be construed as an endorsement or recommendation on the part of the IAEA.*

# Contents

<b>1</b>	<b>Background</b>	<b>2</b>
<b>2</b>	<b>Installations - Virtual Environments and Software Packages</b>	<b>3</b>
2.1	Miniconda3 (conda) and Mamba . . . . .	3
2.2	Git Installation and Repo Cloning . . . . .	4
2.3	Required Libraries with Mamba . . . . .	4
2.3.1	Automatically (faster) . . . . .	5
2.3.2	Manually (slower) . . . . .	5
<b>3</b>	<b>Running Jupyter</b>	<b>6</b>
3.1	Editing the Configuration File . . . . .	6
3.2	Running a RCNV_seq-template Jupyter Notebook . . . . .	9
3.2.1	Installing Required Libraries (optional) . . . . .	9
3.2.2	Loading Required Libraries (mandatory) . . . . .	9
3.2.3	User Input (mandatory) . . . . .	10
3.2.4	CNV Calculations . . . . .	10
3.2.5	Plotting . . . . .	10
3.2.6	Plotting a Zoomed-In Region of One Chromosome . . . . .	11
<b>4</b>	<b>Example Tutorial - CNV Analysis of Sorghum</b>	<b>12</b>
4.1	Data Download . . . . .	12
4.1.1	Linux Terminal . . . . .	13
4.1.2	Web-Browser . . . . .	13
4.2	General Installations . . . . .	13
4.3	Github Repository Cloning . . . . .	13
4.4	Virtual Environment Creation . . . . .	14
4.5	Configuration File Editing . . . . .	14
4.6	Jupyter Notebook Analysis . . . . .	15
<b>5</b>	<b>References</b>	<b>26</b>

---

# 1 Background

Copy number variation (CNV) analysis using CNVseq, R, Jupyter Notebooks, Miniconda3, Mamba, and Git.

All the commands run in a Linux terminal are preceded by the \$ prompt sign. To run a command, copy/past the command without the \$ sign. Those commands run in a Jupyter Notebook are preceded by the *In [ ]:*

---

**Note:** This is not an official IAEA publication but is made available as working material. The material has not undergone an official review by the IAEA. The views expressed do not necessarily reflect those of the International Atomic Energy Agency or its Member States and remain the responsibility of the contributors. The use of particular designations of countries or territories does not imply any judgement by the publisher, the IAEA, as to the legal status of such countries or territories, of their authorities and institutions or of the delimitation of their boundaries. The mention of names of specific companies or products (whether or not indicated as registered) does not imply any intention to infringe proprietary rights, nor should it be construed as an endorsement or recommendation on the part of the IAEA.

---

## 2 Installations - Virtual Environments and Software Packages

Before installing any necessary software, it is recommended to check if the computer is running 32-bit or 64-bit for downloading Miniconda3. Run the following to verify the system:

```
$ uname -m
```

### 2.1 Miniconda3 (conda) and Mamba



Download the Miniconda3, or simply “conda”, installer:

- [Miniconda3 installer for Linux](#)

Run the downloaded installer (for a 64-bit system):

```
$ bash Miniconda3-latest-Linux-x86_64.sh
```

Open a new terminal window for conda to take effect. The word (*base*) should appear in front of the computer name in the terminal window, like so:

A terminal window with a dark background and a menu bar at the top containing 'File', 'Edit', 'View', 'Search', 'Terminal', and 'Help'. The prompt is '(base) anibal@anza-ThinkPad-T580:~\$' with a white cursor. The rest of the terminal is empty.A terminal window with a dark background and a menu bar at the top containing 'File', 'Edit', 'View', 'Search', 'Terminal', and 'Help'. The prompt is '(base) anibal@anza-ThinkPad-T580:~\$' with a white cursor. The rest of the terminal is empty.

Verify the installation and update conda in new terminal window with:

```
$ conda env list  
$ conda update --all  
$ conda upgrade --all
```

Install mamba library/package manager that will be used for installing software dependencies of the tool:

```
$ conda install mamba --yes
```

## 2.2 Git Installation and Repo Cloning

Git is required for cloning locally (downloading a copy to your local computer) the PBGL CNVseq Github repository. Git and Github are used for version control of software. It keeps track of development, releases, and issues of a software project.

Install **git** for cloning the **pbgl-cnvseq** software repository from Github, where the latest version of the tool resides:

```
$ mamba install git --yes
```

After the installation, clone PBGL's CNVseq repository, **pbgl-cnvseq**, to the local computer in any desired directory.

```
$ git clone https://github.com/pbgl/pbgl-cnvseq.git
```

The cloning process will depict the following:

```
Cloning into 'pbgl-cnvseq'...
remote: Enumerating objects: 628, done.
remote: Counting objects: 100% (336/336), done.
remote: Compressing objects: 100% (239/239), done.
remote: Total 628 (delta 109), reused 292 (delta 78), pack-reused 292
Receiving objects: 100% (628/628), 11.79 MiB | 7.03 MiB/s, done.
Resolving deltas: 100% (190/190), done.
```

The **pbgl-cnvseq** repository should have been clones successfully. Verify that the download is complete by listing the folders/files in the directory.

```
$ ls -l
```

The folder called **pbgl-cnvseq** should be listed in the directory.

## 2.3 Required Libraries with Mamba

CNVseq has multiple dependencies, listed below:

- Samtools
- R
  - configr
  - ggplot2
  - BiocManager
  - Bioconductor-GenomicAlignments
  - Bioconductor-GenomeInfoDb
- Jupyter Notebook
  - IRkernel

The necessary R packages are installed through the Jupyter Notebook. It proves as a faster and error-free way to install **configr**, **ggplot2**, **Biocmanager**, **Bioconductor-GenomicAlignments**, and **Bioconductor-GenomeInfoDb** packages.

There are two ways to install the rest of the necessary libraries to run CNV-seq: automatically or manually.

### 2.3.1 Automatically (faster)

One YAML file, **environment.yml**, is provided to automatically create a virtual environment and install the dependent libraries through mamba. The file creates the **cnvseq** virtual environment, along R, Jupyter Notebook, and the R-kernel in Jupyter. It also installs the dependent R libraries. Run **environment.yml**:

```
$ mamba env create --file envs/environment.yml
```

Once done, a list of the virtual environments available can be seen by running:

```
$ conda env list
```

Activate (enter) the recently-created virtual environment **cnvseq**:

```
$ conda activate cnvseq
```

Once done, the virtual environment should be activated and all the necessary packages should be installed. This can be verified with:

```
$ conda list
```

### 2.3.2 Manually (slower)

To manually create and activate an environment, run:

```
$ conda create --name cnvseq
$ conda activate cnvseq
```

Start running the installations of the necessary libraries:

```
$ mamba install --channel conda-forge notebook r-irkernel r-biocmanager --yes
$ mamba install --channel bioconda samtools bioconductor-genomicalignments_
↳bioconductor-genomeinfodb --yes
$ mamba install --channel r r-ggplot2 --yes
$ mamba install --channel pcgr r-configr --yes
```

Once done, all the necessary packages should be installed. This can be verified with:

```
$ conda list
```

## 3 Running Jupyter

To activate Jupyter, run the following in the terminal:

```
$ jupyter notebook
```

This command will start a Jupyter session inside the directory the command is run. The user can navigate between directories, visualize files, and edit files in a web browser by clicking on directories or files, respectively.

Look for the directory **pbgl-cnvsq** and click on it. Click on **tool** directory, which contains three directories and two Jupyter Notebooks. Here is a breakdown of each:

- *config*:
  - directory containing configuration files specifying file paths, parameter definitions, and comparison lists
- *helper-functions*:
  - directory containing R scripts with functions to calculate and plot CNVs
- *output*:
  - directory that will contain both tab-files and images output after running a CNVseq analysis
- two Jupyter Notebooks:
  - RCNV-seq-PBGL-sorghum-analysis-example.ipynb
    - \* example analysis of a comparison between a control and mutant of sorghum
  - RCNV-seq-template.ipynb
    - \* template for the user

---

**Note:** Jupyter lets the user duplicate, rename, move, download, view, or edit files in a web browser. This can be done by clicking the box next to a file and choosing accordingly.

---

### 3.1 Editing the Configuration File

In order to run the CNVseq Jupyter Notebook, the user needs to feed it with a configuration file (**config-CNVseq.yml**) that specifies the paths to the bam files, comparisons to be done, chromosomes to analyze, and parameter definitions for calculating and plotting CNVs.

The configuration file **config-CNVseq.yml** can be found in the **pbgl-cnvsq/tool/config** directory. The configuration file contains the following:

```
# parameters for CNV/window size calculations defaults from HLiag:
# https://github.com/hliang/cnv-seq
parameters:
  annotate: TRUE
  bed_file_present: FALSE
  bigger: 1.5
  log2: 0.6
  pvalue: 0.001
  window_size: 10000

# list of chromosomes to be analyzed
chromosomes:
```

(continues on next page)

```

- Chromosome1
- Chromosome2
- Chromosome3
- AnotherChromosome

# folder to store outputs; should be inside "output" directory, example:
# output_path: output/run-name or output/organism
output_path: output/organism-being-analyzed

# path to bed file for varying window sizes (optional)
bed_path: /path/to/bed/file.bed

# paths to bam files
paths:
  control-1: &control-1 /path/to/bam/file/control-1.bam
  mutant-1: &mutant-1 /path/to/bam/file/mutant-1.bam
  mutant-2: &mutant-2 /path/to/bam/file/mutant-2.bam

# comparisons to be analyzed
comparisons:
  control-1-vs-mutant-1:
    control: *control-1
    mutant: *mutant-1
  control-1-vs-mutant-2:
    control: *control-1
    mutant: *mutant-2

```

---

**Note:** The user needs to edit **config-CNVseq.yml** to point towards bam/bed files; specify comparisons and chromosomes to analyze; and define the parameters to calculate/plot CNVs.

---

One example configuration files is provided (**config-CNVseq-PBGL-sorghum-analysis-example.yml**). The configuration file **config-CNVseq.yml** contains multiple fields to be defined by the user.

- *parameters:*
  - parameters used to create window sizes, thresholds, plots, etc
  - the parameter defaults are provided according to HLiang's original values
- *chromosomes:*
  - list of chromosome names to analyze
  - chromosome names can be found in a bam file's header using the following samtools command:

```
$ samtools view -h my_bam_file.bam | less -S
```

- *output\_path:*
  - directory that will contain both tab-files and images output after running a CNVseq analysis
  - the defined path will be inside the **pbgl-cnvseq/tool/output** directory
  - *images:*
    - \* directory that will store the CNV image outputs per comparison of:
    - \* all chromosomes in one plot
    - \* each chromosome individually in one plot



- \* any zoomed-in region plot of one chromosome
- *tab-files*:
  - \* directory that will contain two types of output tab-delimited files:
    - hits used to calculate CNVs, containing chromosome, start, end, width, control coverage, mutant coverage
    - CNVs per chromosome per comparison, containing CNV number, chromosome, start, end, size, log2, and p-value
- *bed\_path*:
  - path pointing to bed file containing targets if using varying window sizes
  - it is recommended to run the analysis without a **.bed** file; plots will only reflect the targets in a **.bed** file
  - if a **.bed** file is provided, the *bed\_file\_present* parameter under the *parameters* section has to be changed to *TRUE*
- *paths*:
  - sample names and their respective paths to **.bam** files
  - samples can be named as desired but the sample name must be repeated after the colon and prefixed with a & sign
  - the & prefix sign is used to reference the sample's path in different places of the same configuration file
  - example use:

```
paths:
mysample: &mysample /home/username/bam_files/mysample.bam
XYZ-123: &XYZ-123 /home/username/bam_files/XYZ-123.bam
potato95: &potato95 /home/username/bam_files/potato95.bam
```

- *comparisons*:
  - comparison names with respective control and mutant samples per comparison
  - each comparison can be named as desired
  - the sample names to be used as *control* and *mutant* need to be prefixed by a \* sign
  - the \* prefixed sign is used to extract the sample's path defined in the *paths* section
  - example:

```
comparisons:
comparison-1:
  control: *mysample
  mutant: *potato95
a-different-comparison-278asd:
  control: *mysample
  mutant: *XYZ-123
```

## 3.2 Running a RCNV\_seq-template Jupyter Notebook

---

**Note:** It is recommended to duplicate the **RCNV-seq-template** notebook and then renaming the copy before doing any edits to the notebook.

---

In the **pbgl-cvnsseq/tool** directory, click on **RCNV-seq-template** and a new tab in your web-browser will open the notebook.

The notebook contains cells that are populated by text or code. Instructions are provided in the notebook to guide the user. To run a cell, click on the corresponding cell and click on the *Run* button on the top of the notebook. Another way to run a cell can be done by clicking on the corresponding cell and pressing **Ctrl + Enter** or **Shift + Enter**.

The notebook consists of 6 sections:

1. Installing Required Libraries (optional)
2. Loading Required Libraries (mandatory)
3. User Input (mandatory)
4. CNV Calculations
5. Plotting
6. Plotting a Zoomed-In Region of One Chromosome

### 3.2.1 Installing Required Libraries (optional)

- libraries being installed:

```
In [ ]: # install necessary libraries using R functions
        if (!requireNamespace("BiocManager", quietly = TRUE))
            install.packages("BiocManager")

        BiocManager::install(c("GenomicAlignments", "GenomeInfoDb"))
        install.packages("configr")
        install.packages("ggplot2")
```

- to be run if the mamba installations were not successful or the loading of the required libraries fails under the **Loading Required Libraries** section

### 3.2.2 Loading Required Libraries (mandatory)

- libraries being loaded:

```
In [ ]: # load necessary libraries
        library(GenomicAlignments)
        library(ggplot2)
        library(configr)

        # specify source R script with helper functions
        source("helper-functions/RCNV_seq-helper.R")
        source("helper-functions/cnvHLiang.R")
```

- this cell will load necessary libraries and scripts containing the necessary functions to be used
- if running this cell fails, some libraries may be missing from the installation

- to fix this issue, run the installations under the **Installing Required Libraries**

### 3.2.3 User Input (mandatory)

- the user needs to write the appropriate name of the configuration file being used in the following cell:

```
In [ ]: configPath <- "config/config-CNVseq.yml"
```

- the bottom cell will extract the fields defined in the configuration file:

```
In [ ]: config <- read.config(configPath)
```

### 3.2.4 CNV Calculations

- function to calculate all the hits and CNVs:

```
In [ ]: cnvCalculate(config)
```

- output tabulated files are stored inside **pbgl-cnvseq/tool/output/name\_defined\_in\_output\_path\_of\_config/tab-files** directory

### 3.2.5 Plotting

- function to plot two types of images:
  1. CNVs of all chromosomes in the same plot
  2. CNVs of one chromosome per plot

```
In [ ]: cnvPlot(config, imgType="", yMin= , yMax= )
```

- function parameters are:
  - *config* - configuration file defined under **User Input** section
  - *imgType* - image extension to use; available options are: *png*, *jpeg*, *svg*, and *pdf*; default to *png*
  - *yMin* - y-axis bottom limit; default to -5
  - *yMax* - y-axis upper limit; default to 5
- output images are stored inside the **pbgl-cnvseq/tool/output/name\_defined\_in\_output\_path\_of\_config/images** directory
- it is recommended to inspect the CNV-plots with default y-limits and then modify
- *cnvPlotting* can be used in two ways: with default values or user-defined values
  - with default values, the parameters can be omitted and set to *imgType="png"*, *yMin=-5* and *yMax=5*
  - the following two commands have the same parameter values and will output the same plots:

```
In [ ]: cnvPlot(config)
        cnvPlot(config, imgType="png", yMin=-5, yMax=5)
```

### 3.2.6 Plotting a Zoomed-In Region of One Chromosome

- function to plot a specific zoomed-in region of one chromosome

```
In [ ]: cnvPlotZoom(config, tabFile, chromosome="", start= , end= , yMin= , yMax= ,  
↪imgType="")
```

- function parameters are:
  - *config* - configuration file defined under **User Input** section
  - *tabFile* - tabulated file containing all hits; requires user-input to define path to all-hits tabulated file
  - *chromosome* - chromosome name to focus on; default to *NA*
  - *start* - start of window in bp; both scientific notation is accepted: *100000* or *10e4*; defaults to *NA*
  - *end* - end of window in bp; the same applies as in the *start* parameter; defaults to *NA*
  - *yMin* - y-axis bottom limit; default to *-5*
  - *yMax* - y-axis upper limit; default to *5*
  - *imgType* - image extension to use; available options are: *png*, *jpeg*, *svg*, and *pdf*; defaults to *png*
- requires path definition of *tabFile* parameter in the cell:

```
In [ ]: tabFile <- "output/run-name/tab-files/tab-file.tab"
```

## 4 Example Tutorial - CNV Analysis of Sorghum

In this section of the manual, an example analysis of sorghum will be shown in a step-by-step process. The data has been subset to three chromosomes Chr04, Chr05, and Chr09. The CNVseq analysis will depict a deletion on chromosome Chr09.

The tutorial is divided between the following sections:

1. Data Download
2. General Installations
3. Github Repository Cloning
4. Virtual Environment Creation
5. Configuration File Editing
6. Jupyter Notebook Analysis

### 4.1 Data Download

Running this tutorial requires two binary alignment map (BAM) files of sorghum crop: one control and one mutant. Each **.bam** file is subset for chromosomes Chr04, Chr05, and Chr09; this is reflected on their names. They have been both sorted and indexed using samtools.

The following links can be used to download the necessary **.bam** and **.bam.bai** (index) files:

- con-2\_S1 (sorghum control)
  - [https://bss1innov1nafa1poc1.blob.core.windows.net/sample-container/2021\\_Training/con-2\\_S1-Chromes-04-05-09.bam](https://bss1innov1nafa1poc1.blob.core.windows.net/sample-container/2021_Training/con-2_S1-Chromes-04-05-09.bam)
  - [https://bss1innov1nafa1poc1.blob.core.windows.net/sample-container/2021\\_Training/con-2\\_S1-Chromes-04-05-09.bam.bai](https://bss1innov1nafa1poc1.blob.core.windows.net/sample-container/2021_Training/con-2_S1-Chromes-04-05-09.bam.bai)
- D2-1\_S7 (sorghum mutant)
  - [https://bss1innov1nafa1poc1.blob.core.windows.net/sample-container/2021\\_Training/D2-1\\_S7-Chromes-04-05-09.bam](https://bss1innov1nafa1poc1.blob.core.windows.net/sample-container/2021_Training/D2-1_S7-Chromes-04-05-09.bam)
  - [https://bss1innov1nafa1poc1.blob.core.windows.net/sample-container/2021\\_Training/D2-1\\_S7-Chromes-04-05-09.bam.bai](https://bss1innov1nafa1poc1.blob.core.windows.net/sample-container/2021_Training/D2-1_S7-Chromes-04-05-09.bam.bai)

There are two additional ways to download the **.bam** files and their respective indices, besides clicking the links above:

1. Linux terminal
2. Web-Browser

### 4.1.1 Linux Terminal

Open a new terminal and navigate to a directory of choice. We recommend creating a directory to store the data and running *wget* in the respective location:

```
$ mkdir bam_files
$ cd bam_files
$ wget https://bsslinnovlnafalpocl.blob.core.windows.net/sample-container/2021_
↪Training/con-2_S1-Chromes-04-05-09.bam
$ wget https://bsslinnovlnafalpocl.blob.core.windows.net/sample-container/2021_
↪Training/con-2_S1-Chromes-04-05-09.bam.bai
$ wget https://bsslinnovlnafalpocl.blob.core.windows.net/sample-container/2021_
↪Training/D2-1_S7-Chromes-04-05-09.bam
$ wget https://bsslinnovlnafalpocl.blob.core.windows.net/sample-container/2021_
↪Training/D2-1_S7-Chromes-04-05-09.bam.bai
```

### 4.1.2 Web-Browser

Open a web-browser of preference. Copy/paste the links provided above in the address bar. This should automatically begin the download. Move the downloaded files to a location of personal preference.

## 4.2 General Installations

Open a web browser and copy/paste the following link to download Miniconda3:

- <https://docs.conda.io/en/latest/miniconda.html#linux-installers>

After download, open a new terminal window, navigate to the directory with the downloaded Miniconda3 installer, and run the installation.

---

**Note:** The downloaded Miniconda3 installer file might not match the one run in this example. Please, type the corresponding name of the **.sh** file downloaded.

---

```
$ bash Miniconda3-latest-Linux-x86_64.sh
```

Once the Miniconda3 installation is done, close the terminal and open a new one. The word (*base*) should be present to the left of the computer name in the prompt. Update conda and install mamba.

```
$ conda update --all --yes
$ conda upgrade --all --yes
$ conda install mamba --yes
```

### 4.3 Github Repository Cloning

Git or Github is used for storage and version control of software projects. Git is used to manage in a local Linux terminal. First, git will need to be installed. In a new terminal window, run the installation command:

```
$ mamba install git --yes
```

After installing git, create a location to store Github repositories, navigate into it, and clone (download a copy locally) of PBGL's **pbgl-cnvsseq** repository.

```
$ mkdir Github
$ cd Github
$ git clone https://github.com/pbgl/pbgl-cnvsseq
```

A successful cloning will output the following:

```
Cloning into 'pbgl-cnvsseq'...
remote: Enumerating objects: 628, done.
remote: Counting objects: 100% (336/336), done.
remote: Compressing objects: 100% (239/239), done.
remote: Total 628 (delta 109), reused 292 (delta 78), pack-reused 292
Receiving objects: 100% (628/628), 11.79 MiB | 7.03 MiB/s, done.
Resolving deltas: 100% (190/190), done.
```

Navigate into the cloned repository,

```
$ cd pbgl-cnvsseq
```

## 4.4 Virtual Environment Creation

Once inside the **pbgl-cnvsseq** directory, create the **cnvsseq** virtual environment, install the necessary libraries, and activate the newly created **cnvsseq** virtual environment.

```
$ mamba env create --file envs/environment.yml
$ conda activate cnvsseq
```

The name (*cnvsseq*) environment should be reflected to the left of the computer name in the terminal command prompt.

## 4.5 Configuration File Editing

Open a Jupyter session by running,

```
$ jupyter notebook
```

Navigate to **pbgl-cnvsseq/tool/config** and click on **config-CNVseq-PBGL-sorghum-analysis-example.yml**. This will open the configuration file in a new web-browser tab. Copy/paste the following in the configuration file.

---

**Note:** The paths to the bam files will not match. They need to be edited accordingly to point towards the location of the bam files stored locally in the user's computer.

---

```
# parameters for CNV/window size calculations defaults from HLiag:
# https://github.com/hliag/cnv-seq
parameters:
  annotate: TRUE
  bed_file_present: FALSE
  bigger: 1.5
  log2: 0.6
  pvalue: 0.001
  window_size: 10000

# list of chromosomes to be analyzed
chromosomes:
```

(continues on next page)

(continued from previous page)

```
- Chr04
- Chr05
- Chr09

# folder to store outputs; should be inside "output" directory, example:
# output_path: output/run-name or output/organism
output_path: output/PBGL-sorghum-analysis-example

# path to bed file for varying window sizes (optional)
bed_path:

# paths to bam files
paths:
  con-2_S1: &con-2_S1 /home/anibal/bam_files/sorghum/con-2_S1-Chromes-04-05-09.bam
  D2-1_S7: &D2-1_S7 /home/anibal/bam_files/sorghum/D2-1_S7-Chromes-04-05-09.bam

# comparisons to be analyzed
comparisons:
  con-2_s1-vs-D2-1_S7:
    control: *con-2_S1
    mutant: *D2-1_S7
```

Save the file and close the tab.

## 4.6 Jupyter Notebook Analysis

In the open Jupyter session, navigate to the **pbgl-cnvmseq/tool** directory and click on the **RCNV-seq-template.ipynb** Jupyter Notebook. Begin by running the **Loading Required Libraries** section. This can be done by clicking on the cell to be run followed by clicking the *Run* button on the top of the Jupyter Notebook.

```
In [ ]: # load necessary libraries
library(GenomicAlignments)
library(ggplot2)
library(configr)

# specify source R script with helper functions
source("helper-functions/RCNV_seq-helper.R")
source("helper-functions/cnvHLiang.R")
```

If loading the libraries is unsuccessful (specified in the cell output), an attempt to install the missing libraries can be done by running the **Installing Required Libraries** section,

```
In [ ]: # install necessary libraries using R functions
if (!requireNamespace("BiocManager", quietly = TRUE))
  install.packages("BiocManager")

BiocManager::install(c("GenomicAlignments", "GenomeInfoDb"))
install.packages("configr")
install.packages("ggplot2")
```

Run the **User Input** section after correctly editing the name/path of the configuration file:

```
In [ ]: configPath <- "config/config-CNVseq.yml"
```

Load the contents of the configuration file by running the cell,



```
In [ ]: config <- read.config(configPath)
```

Copy number variants (CNVs) are calculated under the **CNV Calculations** section by running the cell,

```
In [ ]: cnvCalculate(config)
```

The output tabulated files of all-hits and CNVs can be inspected in the **pbgl-cnvseq/tool/output/PBGL-sorghum-analysis-example/tab-files** directory. These can be used to import data to other softwares of choice for downstream data analysis. A successful run has the following output right under the input cell,

```
Comparison: con-2_s1-vs-D2-1_S7
Comparing samples:
/home/anibal/bam_files/sorghum/con-2_S1-Chromes-04-05-09.bam
vs.
/home/anibal/bam_files/sorghum/D2-1_S7-Chromes-04-05-09.bam

[1] "chromosome:  Chr04"
[1] "chromosome:  Chr05"
[1] "chromosome:  Chr09"
[1] "cnv_id:  1  of  25"
[1] "cnv_id:  2  of  25"
[1] "cnv_id:  3  of  25"
[1] "cnv_id:  4  of  25"
[1] "cnv_id:  5  of  25"
[1] "cnv_id:  6  of  25"
[1] "cnv_id:  7  of  25"
[1] "cnv_id:  8  of  25"
[1] "cnv_id:  9  of  25"
[1] "cnv_id: 10  of  25"
[1] "cnv_id: 11  of  25"
[1] "cnv_id: 12  of  25"
[1] "cnv_id: 13  of  25"
[1] "cnv_id: 14  of  25"
[1] "cnv_id: 15  of  25"
[1] "cnv_id: 16  of  25"
[1] "cnv_id: 17  of  25"
[1] "cnv_id: 18  of  25"
[1] "cnv_id: 19  of  25"
[1] "cnv_id: 20  of  25"
[1] "cnv_id: 21  of  25"
[1] "cnv_id: 22  of  25"
[1] "cnv_id: 23  of  25"
[1] "cnv_id: 24  of  25"
[1] "cnv_id: 25  of  25"

CNV percentage in genome: 1.6%
CNV nucleotide content: 3140000
CNV count: 25
Mean size: 125600
Median size: 85000
Max Size: 750000
Min Size: 20000

cnv chromosome    start    end      size    log2    p.value
CNVR_1           Chr05    66247501  66267500  20000    1.013742    1.222196e-
↪59
CNVR_2           Chr09    23497501  23527500  30000    -2.288917    3.556107e-
↪238
```

(continues on next page)

(continued from previous page)

CNVR_3	Chr09	23547501	23917500	370000	-2.684822	0
CNVR_4	Chr09	23927501	24132500	205000	-2.639143	0
CNVR_5	Chr09	24142501	24302500	160000	-2.359219	0
CNVR_6	Chr09	24337501	24417500	80000	-1.907744	0
CNVR_7	Chr09	24437501	24467500	30000	-0.912225	1.86383e-
↪68						
CNVR_8	Chr09	24512501	24652500	140000	-2.799081	0
CNVR_9	Chr09	24657501	24692500	35000	-2.466217	9.479478e-
↪297						
CNVR_10	Chr09	24747501	24887500	140000	-2.46515	0
CNVR_11	Chr09	24917501	24942500	25000	-1.980229	9.759448e-
↪172						
CNVR_12	Chr09	24947501	25007500	60000	-2.378318	0
CNVR_13	Chr09	25017501	25102500	85000	-2.464295	0
CNVR_14	Chr09	25107501	25132500	25000	-2.254819	4.908193e-
↪196						
CNVR_15	Chr09	25137501	25202500	65000	-2.90878	0
CNVR_16	Chr09	25212501	25422500	210000	-3.111295	0
CNVR_17	Chr09	25427501	25512500	85000	-2.883919	0
CNVR_18	Chr09	25522501	25637500	115000	-2.617425	0
CNVR_19	Chr09	25647501	25722500	75000	-2.367854	0
CNVR_20	Chr09	25727501	26477500	750000	-2.388962	0
CNVR_21	Chr09	26482501	26507500	25000	-1.727876	9.894624e-
↪147						
CNVR_22	Chr09	26512501	26552500	40000	-1.624917	3.444569e-
↪216						
CNVR_23	Chr09	26562501	26787500	225000	-2.290655	0
CNVR_24	Chr09	26817501	26937500	120000	-1.704837	0
CNVR_25	Chr09	26947501	26972500	25000	-1.303558	1.960969e-
↪100						

The output tabulated files under **pbgl-cnvrseq/tool/output/PBGL-sorghum-analysis-example/tab-files** are the following:

- con-2\_s1-vs-D2-1\_S7-all-chroms-CNVs.tab
- con-2\_s1-vs-D2-1\_S7-window-10000-all-hits.tab

A depiction of the files are shown below,

After the CNV calculations are done, plot the CNVs under the **Plotting** section. Here, the parameters will be set to the following,

```
In [ ]: cnvPlot(cfg, imgType="png", yMin=-5, yMax=5)
```

The output images under **pbgl-cnvrseq/tool/output/PBGL-sorghum-analysis-example/images** are the following:

- con-2\_s1-vs-D2-1\_S7-chromosome-all-window-10000.png
- con-2\_s1-vs-D2-1\_S7-chromosome-Chr04-window-10000.png
- con-2\_s1-vs-D2-1\_S7-chromosome-Chr05-window-10000.png
- con-2\_s1-vs-D2-1\_S7-chromosome-Chr09-window-10000.png

The output plots are the following,

As seen above, each plot contains:

- a title of the chromosome being analyzed

```
1 chromosome start end width test ref
2 Chr04 1 10000 10000 483 521
3 Chr04 5001 15000 10000 590 574
4 Chr04 10001 20000 10000 588 599
5 Chr04 15001 25000 10000 300 342
6 Chr04 20001 30000 10000 241 281
7 Chr04 25001 35000 10000 463 471
8 Chr04 30001 40000 10000 627 627
9 Chr04 35001 45000 10000 507 524
10 Chr04 40001 50000 10000 385 364
11 Chr04 45001 55000 10000 345 290
12 Chr04 50001 60000 10000 381 338
13 Chr04 55001 65000 10000 345 334
14 Chr04 60001 70000 10000 327 338
15 Chr04 65001 75000 10000 345 314
16 Chr04 70001 80000 10000 346 269
17 Chr04 75001 85000 10000 337 292
18 Chr04 80001 90000 10000 362 302
19 Chr04 85001 95000 10000 391 320
20 Chr04 90001 100000 10000 338 320
21 Chr04 95001 105000 10000 294 314
22 Chr04 100001 110000 10000 253 299
23 Chr04 105001 115000 10000 236 283
24 Chr04 110001 120000 10000 314 304
25 Chr04 115001 125000 10000 390 324
26 Chr04 120001 130000 10000 401 314
27 Chr04 125001 135000 10000 404 325
28 Chr04 130001 140000 10000 379 373
29 Chr04 135001 145000 10000 335 352
30 Chr04 140001 150000 10000 335 271
```

Fig. 1: Tabulated file containing all hits of all chromosomes

The screenshot shows a JupyterLab interface with a single tab titled 'con-2\_s1-vs-D2-1\_S7-all-chr'. The browser address bar shows 'localhost:8888/edit/tool/output/PBGL-sorghum-analysis-example/tab-files/co...'. The JupyterLab header includes the Jupyter logo, the file name 'con-2\_s1-vs-D2-1\_S7-all-chr', a checkmark, and the text 'an hour ago'. A 'Logout' button is in the top right. Below the header is a menu bar with 'File', 'Edit', 'View', and 'Language'. On the right of the menu bar is a 'Plain Text' button. The main area contains a code editor with a tabular file of CNVs. The table has 7 columns: 'cnv', 'chromosome', 'start', 'end', 'size', 'log2', and 'p.value'. The rows are numbered 1 to 27. The data is as follows:

cnv	chromosome	start	end	size	log2	p.value
CNVR_1	Chr05	66247501	66267500	20000	1.013742	1.222196e-59
CNVR_2	Chr09	23497501	23527500	30000	-2.288917	3.556107e-238
CNVR_3	Chr09	23547501	23917500	370000	-2.684822	0
CNVR_4	Chr09	23927501	24132500	205000	-2.639143	0
CNVR_5	Chr09	24142501	24302500	160000	-2.359219	0
CNVR_6	Chr09	24337501	24417500	80000	-1.907744	0
CNVR_7	Chr09	24437501	24467500	30000	-0.912225	1.86383e-68
CNVR_8	Chr09	24512501	24652500	140000	-2.799081	0
CNVR_9	Chr09	24657501	24692500	35000	-2.466217	9.479478e-297
CNVR_10	Chr09	24747501	24887500	140000	-2.46515	0
CNVR_11	Chr09	24917501	24942500	25000	-1.980229	9.759448e-172
CNVR_12	Chr09	24947501	25007500	60000	-2.378318	0
CNVR_13	Chr09	25017501	25102500	85000	-2.464295	0
CNVR_14	Chr09	25107501	25132500	25000	-2.254819	4.908193e-196
CNVR_15	Chr09	25137501	25202500	65000	-2.90878	0
CNVR_16	Chr09	25212501	25422500	210000	-3.111295	0
CNVR_17	Chr09	25427501	25512500	85000	-2.883919	0
CNVR_18	Chr09	25522501	25637500	115000	-2.617425	0
CNVR_19	Chr09	25647501	25722500	75000	-2.367854	0
CNVR_20	Chr09	25727501	26477500	750000	-2.388962	0
CNVR_21	Chr09	26482501	26507500	25000	-1.727876	9.894624e-147
CNVR_22	Chr09	26512501	26552500	40000	-1.624917	3.444569e-216
CNVR_23	Chr09	26562501	26787500	225000	-2.290655	0
CNVR_24	Chr09	26817501	26937500	120000	-1.704837	0
CNVR_25	Chr09	26947501	26972500	25000	-1.303558	1.960969e-100

Fig. 2: Tabulated file containing all CNVs of all chromosomes




Fig. 3: CNV plot of all chromosomes




Fig. 4: CNV plot of chromosome Chr04




Fig. 5: CNV plot of chromosome Chr05




Fig. 6: CNV plot of chromosome Chr09



- a sub-title of the comparison being done
- Log2 Ratio in the y-axis
- chromosome name on the x-axis in plot of all chromosomes
- base-pair (bp) position on the x-axis in plot of one chromosome
- a label under the x-axis specifying the location of the stored image

The CNV plot of chromosome Chr09 shows an artifact. Moreover, this artifact seems to lie between base-pairs 2e07 and 3e07. To zoom in, the **Plotting a Zoomed-In Region of One Chromosome** section is provided.

Under the **Plotting a Zoomed-In Region of One Chromosome** section, first specify the location of the all-hits tabulated file of chromosome Chr09,

```
In [ ]: tabFile <- "output/PBGL-sorghum-analysis-example/tab-files/con-2_s1-vs-D2-1_
↳S7-window-10000-all-hits.tab"
```

The run the zooming-in function with the following parameters,

```
In [ ]: cnvPlotZoom(config, tabFile, chromosome="Chr09", start=2.3e7 , end=2.8e7 ,
↳yMin=-5 , yMax=2.5, imgType="png")
```

The output plot, which is also stored in the **pbgl-cnvseq/tool/output/PBGL-sorghum-analysis-example/images**, shows a zoomed-in region of chromosome Chr09 of sorghum. The output path is also specified under the x-axis label.

Knowing where artifacts like this one is located plays a major role in identifying insertions or deletions in an organism. The user can, in turn, use different tools, like the Integrative Genomics Viewer (IGV), for further analysis.

This culminates the tutorial.




Fig. 7: CNV plot of a zoomed-in region of chromosome Chr09

## 5 References

### **BMC Bioinformatics Publication:**

- CNV-seq, a new method to detect copy number variation using high-throughput sequencing

### **GitHub repositories:**

- [hliang/cnv-seq](#)
- [Bioconductor/copy-number-analysis](#)
- [pbgl/pbgl-cnvseq](#)
- [amora197/pbgl-cnvseq](#)